

Ambalal Sarabhai Enterprises Limited

Registered Office : Shantisadan, Mirzapur Road, Ahmedabad-380001. Telephone : +9179-25507671 / 25507073, Fax : +9179-25507483, E-mail : ase@sarabhai.co.in

Ref. No. : Date: 07.04.2023

To,

BSE Limited,

Phiroze Jeejeebhoy Towers, Dalal Street, Fort Mumbai – 400 001 Scrip Code: 500009

Sub: Intimation of New Product(s) Launch by CoSara Diagnostics Private Limited, a Joint Venture of a subsidiary of Ambalal Sarabhai Enterprises Limited ("the Company")

Dear Sir,

Pursuant to Regulation 30 read with Part B of Schedule III of the SEBI (Listing Obligations and Disclosure Requirements) Regulations, 2015, we wish to inform you that **CoSara Diagnostics Private Limited**, a Joint Venture of a subsidiary of the Company and Co-Diagnostics, Inc, Salt Lake City, Utah, USA has launched new products "*RAPID DAIGNOSTICS PRODUCT RANGE*" as per attached.

You are requested to take the same on your record.

Thanking you,


For Ambalal Sarabhai Enterprises Limited

DAMODAR HARGOVINDBHAI SEJPAL Date: 2023.04.07 1630.41 +0530' Damodar H. Sejpal Company Secretary & Compliance Officer

Encl : as above








RAPID DIAGNOSTICS PRODUCT RANGE







SENSITIVE, FAST & AFFORDABLE DIAGNOSTICS

PRODUCT RANGE



SARAGENE[™] HIV (1/2) Card Test

3rd Generation Kit For Differential Detection of HIV 1 & HIV 2

- Coated Ag:gp41 + gp120 for HIV 1 & gp36 for HIV 2
- Specimen Type: Serum/Plasma
- Specimen Volume: 10µl
- Result Interpretation: 10 15 Minutes
- Storage Temperature: 2-40°C

SARAGENE[™] HBsAg (HEPATM) Card Test

Detection of HBV Surface Ag (Australian Ag)

- Specimen Type: Serum/Plasma
- Specimen Volume: 2 3 Drops(Dropper Provided)
- Detection Limit: 0.5ng/mL
- Result Interpretation: 20 Minutes
- Storage Temperature: 2-40°C

SARAGENE[™] HCV Card Test

Detection of Ab Against HCV Antigen

- Specimen Type: Serum/Plasma
- Specimen Volume: 2 3 Drops (Dropper Provided)
- Result Interpretation: 15 20 Minutes
- Storage Temperature: 2-40°C

SARAGENE[™] Syphilis Card Test

Detection of Ab against Treponema pallidum (Syphllis)

- Specimen Type: Serum/ Plasma
- Specimen Volume: 2 3 drops (Dropper Provided)
- Result Interpretation: 15 20 Minutes
- Storage Temperature: 2-40°C









 $\left| \right\rangle$

ONE STEP RAPID TEST

ONE STEP RAPID TEST

COSARA

SARAGENE™ Malaria (Pf-HRP II/Pv-pLDH) Antigen Card Test

Qualitative Detection & Dedifferentiation of Malaria Ag pLDH and HRP II

- Specimen Type: Whole Blood
- Specimen Volume: 5µl
- Detection Limit: 200 Parasites/ µl and confirm to WHO (FIND) Testing Assessment Protocol.
- Result Interpretation: 20 Minutes
- Storage Temperature: 2-40°C

SARAPLEX[™] Dengue Combi Pak (NS1 Ag+IgG/IgM Ab) Card Test

Qualitative Detection of Dengue Early Ag (NS1) & Ab (IgG/IgM)

- Specimen Type: Serum / Plasma
- Specimen Volume: 2 3 drops (Dropper Provided)
- Result Interpretation: 15-20 Minutes
- Storage Temperature: 2-40°C

SARAGENE[™] Dengue Antigen - NS1 Card Test

Qualitative Detection of Dengue Early Ag (NS1)in Human Blood

- Specimen Type: Serum/Plasma
- Specimen Volume: 2 3 drops (Dropper Provided)
- Result Interpretation: 15-20 Minutes
- Storage Temperature: 2-40°C

SARAGENE[™] Dengue Antibody - IgG / IgM Card Test

Qualitative Detection of Dengue Ab (IgG/IgM)

- Specimen Volume: 2 3 drops (Dropper Provided)
- Result Interpretation: 15-20 Minutes
- Storage Temperature: 2-40°C

SARAGENE[™] Typhoid IgG/IgM Card Test


Detection And Differentiation of Anti

- Salmonella Typhi (S. Typhi) IgG and IgM
- Specimen Type: Serum/ Plasma
- Specimen Volume: 2 3 drops (Dropper Provided)
- Result Interpretation: 15 20 Minutes
- Storage Temperature: 2-40°C













Product Code	Product Name	Pack Size
	SARAGENE™	
RPHI01	HIV (1/2) Card Test	50 Test
RPHB01	HBsAg (HEPATM) Card Test	50 Test
RPHC01	HCV Card Test	50 Test
RPSY01	Syphilis Card Test	50 Test
RPMA01	Malaria (Pf-HRP II/Pv-pLDH) Antigen Card Test	50 Test
RPDE01	Dengue Antigen - NS1 Card Test	25 Test
RPTY01	Typhoid IgG/IgM Card Test	25 Test
RPDE02	Dengue Antibody - IgG / IgM Card Test	25 Test
	SARAPLEX™	
RPC001	Dengue Combi Pak (NS1 Ag+IgG/IgM Ab) Card Test	10 Test



Reg. Office :Shantisadan, Mirzapur Road,Ahmedabad - 380 001, India.www.cosara.inTechnic

Factory Address : Sarabhai Campus,
Opp. Ranoli Railway Station,01, India.Ranoli - 391350, Dist. Vadodara, India.General Inquiries : info@cosara.in
Technical Support: support@cosara.in